

WHERE ID-GENE™ DIATOM MOLECULAR INDEX CAN BE APPLIED?

The index has been developed for rivers and streams in Switzerland. However, it is possible to develop its customized versions that can be adapted to other habitats and geographic regions. Such development of customized tests may require an additional step of calibration for local conditions.

HOW TO START?

If you are interested in applying **ID-Gene™ Diatom Molecular Index** to evaluate the biological quality of watercourses in your area, please contact us and we will provide you with protocols and material for samples preservation.

WHAT WE OFFER?

- rapid service (processing time max 2 weeks for 20-50 biofilm samples)
- high quality control
- competitive prices
- final report comprising diatom index value and a list of species for each sample

References

- Swiss Federal Office for the Environment (FOEN/BAFU) *Méthodes d'analyse et d'appréciation des cours d'eau. Diatomées Niveau R / Methoden zur Untersuchung und Beurteilung der Fliessgewässer*
- Visco J, Apothéloz-Perret-Gentil L, Cordonier A et al. (2015) *Environmental Monitoring: Inferring the Diatom Index from Next-Generation Sequencing Data. Environmental Science & Technology, 49, 7597–7605*
- Apothéloz-Perret-Gentil L, Cordonier A, Straub F, Iseli J, Esling P, Pawlowski J. (2017) *Taxonomy-free molecular diatom index for high-throughput eDNA biomonitoring. Mol Ecol Resour. 17: 1231-1242*

ID-Gene ecodiagnosics

Laboratoire
109 ch. du Pont-du-Centenaire,
1228 Plan-les-Ouates

Office
2a ch. des Hutins, 1232 Confignon
Geneva, Switzerland

Tel: +41 076 693 51 69
info@id-gene.com
www.id-gene.com

A new genetic tool for rapid and cost-effective assessment of the biological quality of watercourses

Why diatoms are important for biomonitoring?

Benthic diatoms are widely used as bioindicators in rivers and streams because of their high sensitivity to environmental changes and well-established taxon-specific ecological tolerances and preferences. The diatoms grow rapidly and respond quickly to changes in chemical, physical, or biological factors. Hence, analyzing the composition of their communities provides an easy method to detect environmental changes due to natural or anthropogenic causes.

Various diatom-based biotic indices have been developed all over the world in order to follow the recommendations of public authorities to protect aquatic ecosystems. In Europe, the Member States shall implement measures decided by Water Framework Directive, which recommends using diatoms for the assessment of ecological status of water bodies. In Switzerland, diatoms, along with macro-invertebrates, macrophytes and fishes, form a part of the Modular Stepwise Procedure for the assessment of watercourses published by the Federal Office for the Environment.

What is the ID-Gene™ Diatom Molecular Index?

Based on sequencing of environmental DNA (eDNA) targeting the diatom community, **ID-Gene** has developed a Diatom Molecular Index as a tool to determine the ecological status of rivers and streams. The diatom DNA is specifically amplified from total eDNA extracted from epilithic biofilm samples and sequenced using high-throughput sequencing technology.

Thousands of diatom DNA sequences are obtained per sample. The analysis of these samples using specifically tailored computer algorithms allows to predict with high accuracy a **Diatom Molecular Index**, as well as to provide a list of referenced diatom species present in the sample.

What are the advantages of using the ID-Gene™ Diatom Molecular Index?

Current biodiversity assessment and monitoring are largely based on morphological identification of bioindicator taxa. However, in the case of diatoms, the microscopic identification of tiny and highly variable diatom frustules is time consuming and requires very good taxonomic expertise, which is not always available. These limiting factors contrast with the need of a fast routine

assessment for the management of water quality. To overcome these issues, we propose the **ID-Gene™ Diatom Molecular Index** to assess water quality directly from diatom high-throughput DNA sequence data. This method allows processing a large number of samples over a short period of time.

HOW DOES THE DIATOM MOLECULAR INDEX WORK?

STEP 1

Biofilm DNA extraction



STEP 2

Diatom PCR amplification



STEP 3

High-throughput sequencing



STEP 4

Data analysis



OUTPUT

Index Value



Species list



HOW DOES OUR INDEX PERFORM IN COMPARISON TO CONVENTIONAL MICROSCOPIC APPROACH?

The **ID-Gene™ Diatom Molecular Index** has been calibrated using the Swiss Diatom Index (DI-CH) that is currently available in Switzerland and recognized by the Swiss Federal Office for the Environment as one of the tools to assess the biological quality of Swiss watercourses. The index involves an autoecological value and a weighting factor, which are calculated from the relative abundance and the repartition of each DNA sequence across the entire dataset.

The cross-validation study conducted on more than 500 samples from Swiss rivers and streams show very good congruence between the values of diatom index inferred from morphology-based and eDNA-based taxonomy-free approaches (Fig A). The difference with the index value given by the traditional method is calculated for each sample and for about three quarter of samples (377/508) the ecological status does not differ by more than 0.5 from the original microscopic assessment (Fig B).

